Service Report



WORK ORDER ID: CWKD6975665

WORK ORDER TYPE: Preventative Maintenance

CUSTOMER: UNIVERSITY OF KS MED CTR-CTSU | DEPARTMENT: 146148-UNIVERSITY OF KS MED CTR-CTSU-BIOMED | ACCOUNT NUMBER: 146148

CONTACT: CUSTOMER PO#:

CONTACT PHONE: CUSTOMER REFERENCE:

BILLING ADDRESS: PO Box 3160 Mail Stop 6011 Kansas City, KS 66103 SHIPPING ADDRESS: 4350 Shawnee Mission Parkway Mail Stop 6011 Fairway, KS 66205

AGILITI ASSET ID: 1887931 | SERIAL #: 0006398 | VENDOR SITE ID: | CUSTOMER ASSET ID:

ASSET TYPE: CENTRIFUGES, TABLETOP, LOW-SPEED, REFRIGERATED | MANUFACTURER: EPPENDORF NORTH AMERICA | MODEL: 5702R

SERVICE REQUESTED: PM DUE

SERVICE PROVIDED TYPE: Scheduled Maintenance

WORK COMPLETED: 01-22-2025

LABOR

SERVICE PROVIDERCATEGORYHOURSPERIOD DEFINITIONSERVICE DATEStephanie DyePM Not Covered0.4M-F 8-501/22/25

SERVICE PROVIDED DETAILS

UNIT PASSED ALL REQUIRED TESTING AT THIS TIME

TEST EQUIPMENT

PREFIX UNIT DESCRIPTION CALIBRATION DATE

TA1 0257 TACHOMETER 04/05/24

CHECKLIST

PROCEDURE: ELECTRICAL SAFETY		
QUESTION		ANSWER
The physical integrity of the power cord, attachment plug, and cord-strain relief shall be confirmed by visual inspection		Pass
Ground Resistance (m#)		250
Leakage Current - System Off, Normal Ground, Closed Neutral (µA)		56
Leakage Current - System Off, Open Ground, Closed Neutral (µA)		N/A
Leakage Current - System On, Normal Ground, Closed Neutral (μA) Leakage Current - System On, Open Ground, Closed Neutral (μA) ECG/EKG Lead Leakage (Equipped New Systems/Post Major Repair): System On, Normal Ground less than 100μA ECG/EKG Lead Leakage (Equipped New Systems/Post Major Repair) System On, Open Ground less than 500μA		N/A N/A
		N/A
		PROCEDURE: EPPENDORF NORTH AMERICA 5702R - ANNUAL
QUESTION		ANSWER
Check the centrifuge and the associated rotors.		Pass
Check that the Acceleration time is < 25 s. Check that the Deceleration time is < 25 s. Overspeed cut-off (test carried out without rotor) = Error 7. Check that the Noise level is < 58 dB (A). Check that the Speed is 4,400 1/min ± 44 1/min. Check timer, 5 minutes, tolerance is ± 2 s. Check temperature - Check that Speed can be set at 1,000 rpm to 4,400 rpm. Check temperature - Check that the temperature is 6 °C ± 2 °C, at room temperature 23 °C ± 2 °C. Check imbalance switch: load rotor with G1 on one side = no error.		PASS PASS
		PASS
		4400
		300
		Pass
		6
		Pass
		Check imbalance switch: load rotor with G2 on one side = Error: Imbalance.
Customer Signature:	Date:	
ouotomer orginature.	Dutc.	
Stanfauit	4/00/05	
Stephanie Dye	1/22/25	
Technician Signature:	Date:	