

# Service Report



**WORK ORDER ID: CWKD6982354**

**WORK ORDER TYPE:** Preventative Maintenance

**CUSTOMER:** UNIVERSITY OF KS MED CTR-CTSU | **DEPARTMENT:** 146148-UNIVERSITY OF KS MED CTR-CTSU-BIOMED | **ACCOUNT NUMBER:** 146148

**CONTACT:**

**CUSTOMER PO#:**

**CONTACT PHONE:**

**CUSTOMER REFERENCE:**

**BILLING ADDRESS:** PO Box 3160 Mail Stop 6011 Kansas City, KS 66103

**SHIPPING ADDRESS:** 4350 Shawnee Mission Parkway Mail Stop 6011 Fairway, KS 66205

**AGILITY ASSET ID :** 1887909 | **SERIAL #:** 20258174 | **VENDOR SITE ID:** | **CUSTOMER ASSET ID:**

**ASSET TYPE:** THERMOMETERS, ELECTRONIC | **MANUFACTURER:** WELCH ALLYN INC | **MODEL:** SURETEMP PLUS 690

**SERVICE REQUESTED:** PM DUE

**SERVICE PROVIDED TYPE:** Scheduled Maintenance

**WORK COMPLETED:** 01-22-2025

## LABOR

SERVICE PROVIDER	CATEGORY	HOURS	PERIOD DEFINITION	SERVICE DATE
Stephanie Dye	PM Not Covered	0.3	M-F 8-5	01/22/25

## SERVICE PROVIDED DETAILS

UNIT PASSED ALL REQUIRED TESTING AT THIS TIME

## CHECKLIST

**PROCEDURE:** WELCH ALLYN SURETEMP PLUS 690 - ANNUAL

**QUESTION**

Perform the Power-On Self Test.

**ANSWER**

Pass

Perform the Model 9600 Plus Calibration Testing procedure. (thermometer's temperature reading and the Tester's Set Point Temperature reading are

within  $\pm 0.2^\circ$  F for each Set Point temperature)  
Calibration Testing. Cal-Key Verification - The temperature should read  $97.3 \pm 0.2^\circ\text{F}$ .

Pass  
97.3

---

Customer Signature:

Date:

*Stephanie Dye*

1/22/25

---

Technician Signature:

Date: