

Blood Pressure Monitor Performance Report

Customer:	KU-CTSU Rainbow	Contact:	Nelda Mallett
Street:	3901 Rainbow Blvd (Sudler 2nd Floor)	Phone #	913-574-3007
City/ST/Zip:	Kansas City, KS 66160	P= Pass F=Fail	

MFG	Model#	S/N	Electrical Safety Check	Condition Controls/ Switches	Condition Power Cord & Strin Relief	Hardware Integrity	Indecators & Display Operation	MFG. Labeling Intact	Normal		Hypertonic		Hypotonic		Oral Themp (F)		SPO2	Date Performed
									120/80	P 60	200/150	P130	60/30	P40	98.6	97%		
AXIA Surgical	TRIA	A17044501148	P	P	P	P	P	P	121/82	60	202/152	130	60/31	40	NA	98%	1/25/2023	
AXIA Surgical	TRIA	A17044501115	P	P	P	P	P	P	121/81	60	202/152	130	61/31	40	NA	98%	1/25/2023	
AXIA Surgical	TRIA	A15094700025	P	P	P	P	P	P	121/82	60	202/151	130	60/31	40	NA	98%	1/25/2023	
Welch Allyn	4400	(21)100005910422	P	P	P	P	P	P	121/81	60	201/152	130	61/30	40	98.6	97%	1/25/2023	
Welch Allyn	4400	(21)1000088334021	P	P	P	P	P	P	121/81	60	202/148	130	62/30	40	98.6	98%	1/25/2023	
Welch Allyn	4400	(21)1000088374021	P	P	P	P	P	P	121/81	60	202/149	130	61/32	40	98.8	98%	1/25/2023	
Welch Allyn	4400	(21)100005710422	P	P	P	P	P	P	122/81	60	203/148	130	62/32	40	98.7	97%	1/25/2023	
Welch Allyn	4400	(21)000088734021	P	P	P	P	P	P	122/81	60	202/148	130	62/31	40	98.6	97%	1/25/2023	
Welch Allyn	4400	(21)100034081822	P	P	P	P	P	P	122/81	60	202/148	130	60/31	40	98.6	97%	1/25/2023	
Welch Allyn	297	200819211	P	P	P	P	P	P	122/82	60	201/152	130	61/32	40	NA	98%	1/25/2023	
Welch Allyn	6000	(21)100045103216	P	P	P	P	P	P	121/81	60	203/149	130	62/32	40	NA	NA	1/25/2023	
DELP G040																		
Welch Allyn	42NTB	201907174	P	P	P	P	P	P	121/81	60	202/151	130	60/32	40	98.1	98%	1/25/2023	
AXIA Surgical	TRIA	A18044501468	P	P	P	P	P	P	121/82	60	202/151	130	60/31	40	NA	98%	1/25/2023	

Technicians	JC	Technician: Jim Cowart	Test equipment used FLUKE- ProSim 8 S#399003 Expires:09/07/2023
Initials		Date: 1/25/2023	
Notes:		Next Due Date: 01/2024	
Calibration Performed by Cobb Refrigeration and Lab Service 8071 Flint St. Lenexa KS, 66212 913-403-0088			

Blood Presure Monitor Performance Report

Customer:		KU-CTSU Rainbow			Contact:		Nelda Mallett										
Street:		3901 Rainbow Blvd (Sudler 2nd Floor)			Phone #		913-574-3007										
City/ST/Zip:		Kansas City, KS 66160					P= Pass F=Fail										
MFG	Model#	S/N	Electrical Safety Check	Condition Controls/ Switches	Condition Power Cord & Strin Relief	Hardware Integrity	Indecators & Display Operation	MFG. Labeling Intact	Normal		Hypertonic		Hypotonic		Oral Themp (F)	SPO2	Date Performed
									120/80	P 60	200/150	P130	60/30	P40	98.6	97%	

Blood Pressure Monitor Performance Report

Customer:		KU-CTSU Rainbow			Contact:		Nelda Mallett										
Street:		3901 Rainbow Blvd (Sudler 2nd Floor)			Phone #		913-574-3007										
City/ST/Zip:		Kansas City, KS 66160					P= Pass F=Fail										
MFG	Model#	S/N	Electrical Safety Check	Condition Controls/ Switches	Condition Power Cord & Strin Relief	Hardware Integrity	Indecators & Display Operation	MFG. Labeling Intact	Normal		Hypertonic		Hypotonic		Oral Themp (F)	SPO2	Date Performed
									120/80	P 60	200/150	P130	60/30	P40	98.6	97%	