

Blood Presure Monitor Performance Report

Customer:	KU-CTSU Rainbow	Contact:	Shelby Williams
Street:	3901 Rainbow Blvd (Sudler 2nd Floor)	Phone #	913-574-3018
City/ST/Zip:	Kansas City, KS 66160	P= Pass F=Fail	

MFG	Model#	S/N	Electrical Safety Check	Condition Controls/ Switches	Condition Power Cord & Strin Relief	Hardware Integrity	Indecators & Display Operation	MFG. Labeling Intact	Normal		Hypertonic		Hypotonic		SPO2	Date Performed
									120/80	P 60	200/150	P130	60/30	P40		
AXIA Surgical	TRIA	A18014600076	P	P	P	P	P	P								
AXIA Surgical	TRIA	A18044501442	P	P	P	P	P	P	120/81	60	202/149	130	59/31	40	98%	1/31/2022
AXIA Surgical	TRIA	A17044501148	P	P	P	P	P	P	120/82	60	202/153	130	61/30	40	97%	1/31/2022
AXIA Surgical	TRIA	A17044501115	P	P	P	P	P	P	120/81	60	202/151	130	59/31	40	98%	1/31/2022
AXIA Surgical	TRIA	A17044501158	P	P	P	P	P	P	121//81	60	203/150	130	60/32	40	98%	1/31/2022
AXIA Surgical	TRIA	A15094700025	P	P	P	P	P	P	121/81	60	204/151	130	61/30	40	98%	1/31/2022
AXIA Surgical	TRIA	A15094700021	P	P	P	P	P	P	121/82	NA	201/152	NA	60/32	NA	NA	1/31/2022
AXIA Surgical	TRIA	A17044501067	P	P	P	P	P	P	120/81	60	202/151	130	60/32	41	97%	1/31/2022
AXIA Surgical	TRIA	A19014800070	P	P	P	P	P	P	119/82	60	201/153	130	61/30	40	97%	1/31/2022
Welch Allyn	297	200819211	P	P	P	P	P	P	121/81	60	202/150	130	61/30	40	97%	1/31/2022
DELP G040																
Welch Allyn	42NTB	201907174	P	P	P	P	P	P	121/81	60	201/151	130	61/32	40	97%	1/31/2022
AXIA Surgical	TRIA	A18044501468	P	P	P	P	P	P	119/81	60	201/151	130	59/31	40	97%	1/31/2022

Technicians JC Technician: Jim Cowart Test equipment used FLUKE- ProSim 8 S#399003 Expires:09/10/2022

Blood Pressure Monitor Performance Report

Customer: KU-CTSU Rainbow		Contact: Shelby Willams															
Street: 3901 Rainbow Blvd (Sudler 2nd Floor)		Phone # 913-574-3018															
City/ST/Zip: Kansas City, KS 66160		P= Pass F=Fail															
MFG	Model#	S/N	Electrical Safety Check	Condition Controls/ Switches	Condition Power Cord & Strin Relief	Hardware Integrity	Indecators & Display Operation	MFG. Labeling Intact	Normal		Hypertonic		Hypotonic		SPO2	Date Performed	
									120/80	P 60	200/150	P130	60/30	P40	97%		
Initials				Date: 01/31/2022													

Notes: **Next Due Date: 01/2023**

Calibration Performed by Cobb Refrigeration and Lab Service 8071 Flint St. Lenexa KS, 66212 913-403-0088

Blood Pressure Monitor Performance Report

Customer: KU-CTSU Rainbow		Contact: Shelby Williams														
Street: 3901 Rainbow Blvd (Sudler 2nd Floor)		Phone # 913-574-3018														
City/ST/Zip: Kansas City, KS 66160		P= Pass F=Fail														
MFG	Model#	S/N	Electrical Safety Check	Condition Controls/ Switches	Condition Power Cord & Strin Relief	Hardware Integrity	Indecators & Display Operation	MFG. Labeling Intact	Normal		Hypertonic		Hypotonic		SPO2	Date Performed
									120/80	P 60	200/150	P130	60/30	P40	97%	

Blood Pressure Monitor Performance Report

Customer: KU-CTSU Rainbow		Contact: Shelby Williams														
Street: 3901 Rainbow Blvd (Sudler 2nd Floor)		Phone # 913-574-3018														
City/ST/Zip: Kansas City, KS 66160		P= Pass F=Fail														
MFG	Model#	S/N	Electrical Safety Check	Condition Controls/ Switches	Condition Power Cord & Strin Relief	Hardware Integrity	Indecators & Display Operation	MFG. Labeling Intact	Normal		Hypertonic		Hypotonic		SPO2	Date Performed
									120/80	P 60	200/150	P130	60/30	P40	97%	