

Service Report



WORK ORDER ID: CWKD4070817

WORK ORDER TYPE: Preventative Maintenance

CUSTOMER: UNIV OF KS MED CTR-CTSU RNBW | **DEPARTMENT:** 177820-UNIV OF KS MED CTR-CTSU RNBW-ROUNDS CTSU RNBW | **ACCOUNT NUMBER:** 177820

CONTACT:

CUSTOMER PO#:

CONTACT PHONE:

CUSTOMER REFERENCE:

BILLING ADDRESS: 3901 Rainbow Blvd MS 1077 2012 Suder Kansas City, KS 66160 | **SHIPPING ADDRESS:** 3901 Rainbow Blvd MS 1077 2012 Suder Kansas City, KS 66160

AGILITY ASSET ID : 1887873 | **SERIAL #:** 100005710422 | **VENDOR SITE ID:** | **CUSTOMER ASSET ID:**

ASSET TYPE: MONITORS, PHYSIOLOGIC, VITAL SIGNS | **MANUFACTURER:** WELCH ALLYN INC A SUBSIDIARY OF HILLROM HOLDINGS INC | **MODEL:** SPOT VITAL SIGNS 4400

SERVICE REQUESTED: Asset Onboarding - Initial Preventative Maintenance for 1887873

SERVICE PROVIDED TYPE: Scheduled Maintenance

WORK COMPLETED: 01-19-2024

LABOR

| SERVICE PROVIDER | CATEGORY | HOURS | PERIOD DEFINITION | SERVICE DATE |
|-------------------|----------------|-------|-------------------|--------------|
| Stephanie Harper | PM Not Covered | 0.22 | M-F 8-5 | 01/19/24 |
| Marcus Boatwright | PM Not Covered | 0.22 | M-F 8-5 | 01/19/24 |

SERVICE PROVIDED DETAILS

Completed SM passed all inspection testing

TEST EQUIPMENT

| PREFIX | UNIT | DESCRIPTION | CALIBRATION DATE |
|--------|------|--------------------|------------------|
| PR3 | 473 | OXIMETER SIMULATOR | 04/14/23 |

CHECKLIST

PROCEDURE: MONITORS, BEDSIDE, PHYSIOLOGIC

QUESTION

ANSWER

| | |
|--|------|
| 01. Check AC plug/line cord/strain relief | Pass |
| 02. Measure chassis ground resistance (milliOhms) | 1 |
| 03. Measure chassis leakage current (microAmps) | 85 |
| 04. Patient Leakage (microAmps) | N/A |
| 05. Check outercase, line cords, rolling stands, wall mounts, modular accessories and interconnecting cables | Pass |
| 06. Check patient interface connections (ECG, IBP, SPO2, Temp, CO2 and NIBP) | Pass |
| 07. Check display panel, feet, screws and handle | Pass |
| 08. Keypad/Control Knob Test | Pass |
| 09. Recorder Test | NA |
| 10. Display Tests | Pass |
| 11. Pixel Test | Pass |
| 12. Color Test | Pass |
| 13. NIBP Tests - Overpressure Voltage Test | Pass |
| 14. NIBP Tests - Static Pressure Calibration | Pass |
| 15. NIBP Tests - Motor Pump Test | Pass |
| 16. NIBP Tests - Leak Test | Pass |
| 17. NIBP Tests - Dump Test (Adult, Neonate) | Pass |
| 18. NIBP Tests - Bleed Rate Test | Pass |
| 19. NIBP Tests - Overpressure Test (Adult, Pediatric, Neonate) | Pass |
| 20. Voltage and Battery Operation Verification | Pass |
| 21. ECG Tests - Initialization | Pass |
| 22. ECG Tests - Leads Off | Pass |
| 23. ECG Tests - Pacer Detect | Pass |
| 24. ECG Tests - Heart Rate | Pass |
| 25. ECG Tests - Alarms | Pass |
| 26. IBP 1 and IBP 2 Verification | NA |
| 27. Temperature Verification | Pass |
| 28. SPO2 Verification | Pass |
| 29. SPO2 Alarm | Pass |
| 30. NIBP Verification | Pass |
| 31. CO2 Operation Verification | NA |

Customer Signature:

Date:

Technician Signature:

Date: