

Service Report



WORK ORDER ID: CWKD8232840

WORK ORDER TYPE: Preventative Maintenance

CUSTOMER: UNIVERSITY OF KS MED CTR-CTSU | **DEPARTMENT:** 146148-UNIVERSITY OF KS MED CTR-CTSU-BIOMED | **ACCOUNT NUMBER:** 146148

CONTACT:

CONTACT PHONE:

BILLING ADDRESS: 4350 Shawnee Mission Parkway Mail Stop 6011 Fairway, KS 66205

CUSTOMER PO#:

CUSTOMER REFERENCE:

SHIPPING ADDRESS: 4350 Shawnee Mission Parkway Mail Stop 6011 Fairway, KS 66205

AGILITI ASSET ID : 1705470 | **SERIAL #:** 20258156 | **VENDOR SITE ID:** | **CUSTOMER ASSET ID:**

ASSET TYPE: OPHTHALMOSCOPE/OTOSCOPES | **MANUFACTURER:** WELCH ALLYN INC | **MODEL:** SURE TEMP 4

SERVICE REQUESTED: Asset Onboarding - Initial Preventative Maintenance for 1705470

SERVICE PROVIDED TYPE: Scheduled Maintenance

WORK COMPLETED: 04-28-2025

LABOR

SERVICE PROVIDER	CATEGORY	HOURS	PERIOD DEFINITION	SERVICE DATE
Stephanie Dye	PM Not Covered	0.4	M-F 8-5	04/28/25

SERVICE PROVIDED DETAILS

UNIT PASSED ALL REQUIRED TESTING AT THIS TIME

CHECKLIST

PROCEDURE: WELCH ALLYN INC A SUBSIDIARY OF HILLROM HOLDINGS INC SURE TEMP 4 MONTHLY

QUESTION

Visually inspect the thermometer for physical damage which might cause future product failure.

ANSWER

Pass

Perform power up display test.	Pass
Perform startup display test.	Pass
Calibration - Verify the values on the display are within the range. (°F)	97.1
Model 9600 Calibration - Verify that the temperatures are within ±0.1° C (±0.2° F).	PASS

<hr/>	<hr/>
Customer Signature:	Date:

<hr/>	<hr/>
Technician Signature:	Date: